

FLAMMABLE

Flammable Material Storage Laboratory Refrigerators and Freezers

- Meet NFPA and OSHA guidelines for flammable material storage in general laboratory work areas
- UL Listed
- Temperature Performance:
 - Refrigerators: 0°C to 10°C, Preset at 4°C
 - Freezers: -18°C to -25°C, Preset at -20°C
 - Dual Temp: Refrigerator Preset at 4°C
Freezer Preset at -20°C

- Exterior Adjustable Mechanical Thermostat
- 120 VAC Operation with Power Cord
- Leveling Legs
- Keyed Door Lock
- Stainless Steel Handle
- White Baked on Powder Coat Finish
- Heavy Duty Powder Coated Shelves
- Durable ABS Material Flat Door Liner

Accessories

- Caster Platform
- Additional Full Shelf



Flammable Series

ABS Model	Cu. Ft	Defrost	Temp	EXT. Doors	INT. Doors	Shelves	Width (Exterior)	Depth	Height	HP	Amps	BTUs	Weight	Shipping Weight
ABT-FR-3004	30	Manual	+4C	1 Solid	0	5 Adjustable	35	35	72	1/3	7	1800	335	359
ABT-FF-3020	30	Manual	-20C	1 Solid	0	3 Fixed/2 Adj	35	35	72	1/3	7	1800	332	356
ABT-FC-30	30	Manual	+4C/-20C	1 Solid	1	3 Adj Refrig	35	35	72	1/3	7	1800	347	374